



## Instructions for Use

**REF** Premium Head & Chin Straps 800-0369

### INTENDED USE

Intended use is to hold patient's head during operative procedure. The intended users of this device are medical professionals within hospitals and surgery centers.

### INSTRUCTIONS

Become familiar with features of patient positioning device before use with patient. Always practice use on a nurse, physician or appropriate volunteer prior to using clinically.



1. Place patient on surgical table/shoulder chair in supine position
2. Once patient is intubated, hold patient's head/neck and raise shoulder chair up so patient is in Fowler (sitting) position
3. Loosen ball joint - adjust head positioner around patient's head making sure patient's ears are uncovered - tighten ball joint
4. Hold head positioner flaps alongside of patient's head, placing Forehead Strap across forehead as shown. Attach securely and avoid covering patient's eyebrows.
5. Center Chin Strap over chin, secure strap to one side of head positioner. Making sure head is properly positioned, secure strap's other end to opposite side of head positioner.



**NOTE:** Head should be securely positioned. If not, make straps a little tighter.

### COMPONENT OVERVIEW

Premium Head and Chin Straps holds a patient's head securely during operative procedure.

### GENERAL SPECIFICATIONS

*Device Dimensions (maximum)*

- Length: 21" +/- 0.5" (53 cm +/- 1 cm) (Head)
- Length: 22" +/- 0.5" (56 cm +/- 1 cm) (Chin)
- Width: 2" +/- 0.5 (5 cm +/- 1 cm) (Both Head and Chin)
- Device Weight: .066 lbs. +/- .02 lbs. (.03 kg +/- .009 kg)
- Attaches to head section of SchureMed Shoulder Chairs
- Single-person installation

### GENERAL INFORMATION

- Product not made with Natural Rubber Latex

- If damaged in shipping, call Customer Service at (888) 724-8763 or (781) 982-7000 to obtain Return Material Authorization number. For product warranty issues, contact Customer Service.
- CE marked medical device according to MDR (EU) 2017/745
- Store device between -4°F to +86°F (-20°C to 30°C)



**SINGLE PATIENT USE—Please dispose after each use**

## USER NOTICE

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## eIFU Language Versions

To download and print the Instructions for Use, please go to [www.schuremed.com](http://www.schuremed.com)

## Symbol Glossary

Symbol	Title
	Manufacturer
	Date of manufacture
	Authorized Representative in the European Community
	Authorized Representative in the Swiss
	Importer
	Serial Number
	Warning
	Medical Device
	Unique Device Identifier
	CE Marking
	Single Patient Use



**SchureMed** (081001460)

452 Randolph Street, Abington, MA 02351 USA

**Toll Free** (888) 724-8763 | **Ph** (781) 982-7000 | **Fax** (781) 982-7001 | [orders@schuremed.com](mailto:orders@schuremed.com)

