



Powered Beach Chair

The only powered surgical beach chair on the market with a self-contained motor

Our double-ball joint head positioning system safely secures the patient's head; even your kyphosis patients and allows patient's head to tilt 30° in every direction. Motion-control technology eliminates manual patient positioning with a simple push-button hand control.

Features

- Keeps the patient and O.R. staff safer and more secure
- Attaches to ANY operating room table with our integrated, adjustable mounting clamps
- 5-Piece Deluxe Pad Set and Lateral Brace included
- 1,350 lbs. of lift securely supports arthroscopy patients up to 500 lbs.
- Shoulder pull-away panels provide excellent surgical exposure
- Articulation range 0-90° in 8 seconds in a slow, controlled movement
- 28"W x 30"D x 47"H, 32 lbs

#800-0004 Powered Beach Chair
#800-0004-EU Powered Beach Chair, EU

Schure Socket XL
#800-0134 US
#800-0134-EU EU
#800-0134-JPN JPN
#800-0134-UK UK
#800-0134-DEN DEN

Optional Accessories
#800-0074 Shoulder Chair Dolly
#800-0050 Multi-Axis Arm Positioner
#800-0164 Disposable Face Masks, 12/cs
#800-0039 Disposable Head & Chin Straps, 12/cs
#508-1345 Reusable Head & Chin Straps, set
#800-0091 Lateral Brace



Patented

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351