



Schure™ Socket XL

Quite possibly the best operating room table clamp you have ever used!

The revolutionary Schure™ Socket XL fits all standard operating room table side rails ($\frac{3}{8}$ "W x $1\frac{1}{8}$ "H). Quickly and easily connects and disconnects securely anywhere along the side rail, even over side rail stand off.

Features

- Hardened clamping surface creates maximum clamping force stopping movement of all accessories
- Accepts any round mounting posts up to $\frac{5}{8}$ " in diameter and rectangular mounting posts up to $1" \times \frac{3}{8}"$
- 800 lb (363 kg) patient weight capacity

Schure™ Socket XL

#800-0134 US
#800-0134-UK UK
#800-0134-EU EU
#800-0134-DEN DEN
#800-0134-JPN JPN

Each clamp sold separately.

*Simplify
your
clamp
inventory*



Patented

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351