



# Schure™ Spine Frame

*A secure, comfortable and adjustable spinal platform*

The Schure™ Spine Frame provides a stable and convenient platform for spine surgeries such as laminectomies, decompressions, and disc surgery. Tapered adjustable pads provide optimal patient safety and comfort.

## Features

- *Adjusts intra-operatively to the desired lordosis to get optimal surgical site exposure*
- *500 lb patient weight capacity*
- *Can be used with any O.R. table and fits the Jackson Spine Table and Allen Flex Frame as well*
- *Frame comes complete with transport storage dolly, crank handle, disposable head rest, and 2-arm cradles*
- *Radiolucent and allows for unrestricted C-arm access*
- *Deluxe foam pads included*

**#800-0235** Schure™ Spine Frame,  
18"W x 7"D x 30"H, 18 lbs

## Optional Accessories


**#508-0373** Disposable Laminectomy Arm Cradles,  
12 prs/cs

**#508-1340** Disposable Richard Slotted Head Rests,  
Extended, 6/cs

**#508-0472** Schure™ Spine Frame Pad Set

**#800-0277** Schure™ Spine Frame Crank Handle

**#800-0074-SF** Dolly, Schure™ Spine Frame



Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

452 Randolph Street, Abington, MA 02351