



Great White® Platinum Stirrups

*Gynecology, laparoscopy
& urology procedures
are safer than ever*

The Great White® Platinum Stirrups provide safe and quick positioning assistance for improved surgical site access while offering exceptionally precise placement.

Features

- Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -30° low to 80° high; Adduction -9° and abduction +25°
- The raised lateral fin boots effectively minimize the possibility of nerve damage
- Included Platinum Boot Pads wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning
- Lithotomy indicators help to accurately reposition legs

#800-0342-PL Great White® Platinum Stirrups

Required Accessory: Purchase Two SpringLoc Clamp

#800-0338 US #800-0338-DEN DEN #800-0338-EU
#800-0338-JPN JPN #800-0338-SWISS SWISS

Optional Accessories

#508-1415 Replacement Stirrup Boots Pads w/Fins
#800-0074-S Great White® Stirrup Dolly
#800-0284 Stirrup Wall Rack
#800-0160 Clear View Leg Drapes



600 lb. (272kg)
Patient Weight
Capacity



Patent
Pending

*Increase your
surgeons' O.R.
performance*

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351