

# Great White® Platinum Stirrups

*Gynecology, laparoscopy  
& urology procedures  
are safer than ever*

The Great White® Platinum Stirrups provide safe & quick positioning assistance for improved surgical site access while offering exceptionally precise placement.

## Features

- Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -25° low to 75° high; Adduction -9° and abduction +25°
- The raised lateral fin boots effectively minimize the possibility of nerve damage
- Included Platinum Boot Pads wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning
- Lithotomy indicators help to accurately reposition legs

#800-0342-PL Great White® Platinum Stirrups

## Required Accessory: Purchase Two SpringLoc Clamp

#800-0338 US  
 #800-0338-EU EU  
 #800-0338-JPN JPN  
 #800-0338-SWISS SWISS



**600 lb. (272 kg)  
Patient Weight  
Capacity**



Patent Pending

## Optional Accessories

#508-1415 Replacement Boots Pads w/Fins  
 #800-0074-S Great White® Stirrup Dolly  
 #800-0284 Stirrup Wall Rack  
 #800-0160 Disp. Clear View Leg Drapes, Sterile, 10/cs\*  
 \*Sterile products are non-returnable.

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351