



## TKR Positioner for Total Knee Replacement

*Incredibly innovative knee positioner that offers a completely secure, dynamic-range of patient positions for TKR procedures from SchureMed*

This is a huge leap forward for the surgeon and patient alike for total knee replacements. The ball-track system at the heel of the leg holder places the operative leg in flexion, extension, tilt and rotation.

### Features

- *Track Locking Bracket securely locks the Base Plate to the table with two thumb screws*
- *Locks in any position with a simple turn of the wrist*
- *Smooth sliding motion for effortless leg positioning*
- *Well leg lies flat – eliminating the need for well-leg positioners, sandbags and gel pads*
- *Easy-to-use and completely rigid once locked*
- *All components are autoclavable*

### TKR Positioner for Total Knee Replacement

**#800-0141** w/24" Track System, 11"W x 8"D x 24"H, 5 lbs

**#800-0171** w/30" Track System, 11"W x 8"D x 30"H, 6 lbs

### CamLoc Clamp

**#800-0271** US      **#800-00271-UK** UK

**#800-0271-EU** EU      **#800-00271-DEN** DEN

### Optional Accessory

**#800-0146** TKR Disp. Sterile TKR Boot Liners, 12/cs



*Increase your surgeons' O.R. performance*

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351