



## TKR Positioner for Total Knee Replacement

*Incredibly innovative knee positioner that offers a completely secure, dynamic-range of patient positions for TKR procedures from SchureMed*

This is a huge leap forward for the surgeon and patient alike for total knee replacements. The ball-track system at the heel of the leg holder places the operative leg in flexion, extension, tilt and rotation.

### Features

- Track Locking Bracket securely locks the Base Plate to the table with two thumb screws
- Locks in any position with a simple turn of the wrist
- Smooth sliding motion for effortless leg positioning
- Well leg lies flat – eliminating the need for well-leg positioners, sandbags and gel pads
- Easy-to-use and completely rigid once locked
- All components are autoclavable

### TKR Positioner for Total Knee Replacement

#800-0141 with 24" L Track System

#800-0171 with 30" L Track System

### CamLoc Clamp

#800-0271 US

#800-00271-UK UK

#800-0271-EU EU

#800-00271-DEN DEN

#800-0271-JPN JPN


#800-00271-SWISS SWISS

### Optional Accessory

#800-0146 TKR Disp. Sterile TKR Boot Liners, 12/cs\*

\*Sterile products are non-returnable.

New Improved Design!



Increase your surgeons' O.R. performance

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351