

Candy Cane Stirrups

Position patients for gynecologic, urologic or laparoscopic procedures effectively

Candy Cane Stirrups are a cost-effective way to position patients for a wide-range of surgical procedures. Lock the medical stirrup's height or change the angle of the patient's leg by simply twisting the ergonomic handles.

Features

- Adjust vertically from 28"-44" (71cm-112cm) with an 8" (20cm) leg clearance
- 18"L x 1.5"W (46cm x 4cm) double-looped, cotton ankle strap
- Lock stirrup height by simply twisting the ergonomic handle
- Connect stirrups to OR table side rail with Schure Socket XL (sold separately)
- Candy Cane Stirrup Pole Pad, Strap Assembly and Hooks included



#800-0012 Candy Cane Stirrups

Required Accessory: Order Two Schure Socket XL

#800-0134 US

#800-0134-EU EU

#800-0134-JPN JPN

#800-0134-UK UK

#800-0134-DEN DEN

Optional Accessories

#790-0006 Candy Cane Strap Replacement, ea

#508-0453 Gel Candy Cane Stirrup Sheath, ea

#512-0003 Candy Cane Repl. Hook, ea

Phone 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351