



Adjustable Height Armboard

Eliminates need for multiple size pads

Accommodates a wide-range of patients and their procedural needs, as well as adjusting vertically to allow the use of x-ray tops without sacrificing a level surface on the O.R. table.

Features

- To adjust armboard vertically, simply place armboard on side rail, loosen handle and raise to correct table height and tighten handle
- Simply pull finger-trigger at end of armboard and quickly move to required position. Release trigger to lock tight.
- Pad sold separately
- 6"W x 26.5"L



Adjustable Height Armboard

#800-0009 US #800-0009-DEN DEN
#800-0009-EU EU #800-0009-JPN JPN
#800-0009-UK UK

Required Accessory: Choose One

#508-0071 2" (5cm) Deluxe Pad
#508-0151 2" (5cm) Softcare Pad

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351