



## Deluxe Rail Clamp

*Accessories can be mounted vertically or horizontally*

The Deluxe Rail Clamp locks all patient positioning accessories and surgical equipment anywhere along the operating table side rail.

### **Features**

- Attaches at OR table rail ends or notches
- Securely locks all accessories with 1" x 3/8" (2.5cm x .95cm) flat mounting post and round mounting posts up to 5/8"D (1.6cm) mounting blades anywhere along the table side rail
- 500 lb (227 kg) patient weight capacity

### **Deluxe Rail Clamp**


**#800-0248** US

**#800-0248-UK** UK

**#800-0248-EU** EU

**#800-0248-DEN** DEN

*Each clamp sold separately.*



**Phone** 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351