



# Great White® Premium Stirrups

*There's a new stirrup in town!*

High quality, feature-rich, and priced right, the Great White® Premium Stirrups are exactly what your surgical team requires

**400 lb. (181 kg)  
Patient Weight  
Capacity**

## Features

- Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -30° low to 80° high; Adduction -9° and abduction +25°
- Included Premium Boot Pads wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning
- Lithotomy indicators help to accurately reposition legs

**#800-0342-PR** Great White® Premium Stirrups w/o Fin

**#800-0342-PRF** Great White® Premium Stirrups w/Fin

## Required Accessory: Purchase Two

Spring Loc Clamp

**#800-0338** US

**#800-0338-EU** EU

**#800-0338-DEN** DEN

**#800-0338-JPN** JPN

**#800-0338-SWISS** SWISS

## Optional Accessories

**#508-1354** Replacement Stirrup Boot Pads w/o Fins

**#508-1415** Replacement Stirrup Boot Pads w/Fins

**#800-0074-S** Dolly, Great White® Stirrup

**#800-0284** Stirrup Wall Rack

**#800-0160** Clear View Leg Drapes



*Increase your  
surgeons' O.R.  
performance*

**Patent  
Pending**

**Phone** 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

452 Randolph Street, Abington, MA 02351

©2018 Schuerch Corporation #065-0064 REV 2.03