



NEW!

Great White® HD 1000 Stirrup Boot System

Discover innovation in our interchangeable stirrup boot system with 1,000 lb Patient Wt. Capacity

Includes 3 Different Style Boots!

Features

- **Proprietary 3-in-1 System** includes Platinum, Bariatric Maxima & Premium Stirrup Boots
- **Stronger!** A new, engineered stainless steel tube increases load capacity to 1,000 lb
- **Redesigned!** 3-axis control handle mimics the human hip's motion for smooth and accurate high/low lithotomy, adduction & abduction positioning
- **Locking Gas Spring** locks stirrup in place—stays locked until trigger is pulled
- **Ranges** Lithotomy Range 25° low to 75° high; Adduction -9° and abduction +25°
- **Included** Premier boot liners wrap the foot, ankle & calf providing safe, comfortable & secure patient positioning
- **Reduced Setup** Quick, easy connection reduces setup time, training and downtime
- **Easily Repositions** Lithotomy indicators accurately reposition legs



Length & lithotomy indicators



Select from the included (A) Bariatric Maxima, (B) Platinum & (C) Premium Boots.



Quick and easy connection reduces setup time making it a cost effective product.

#800-0364 Great White® HD 1000 Stirrup Boot System

Required Accessory: Purchase Two SpringLoc Clamps

#800-0338 US #800-0338-EU EU
 #800-0338-JPN JPN #800-0338-SWISS SWISS

Optional Accessories

- #508-1415 Replacement Boot Pads w/Fins
 - #508-1354 Replacement Boot Pads w/o Fins
 - #508-0559 Replacement Bariatric Maxima Boot Pads
 - #800-0160 Disp. Clear View Leg Drapes, Sterile, 10/cs*
- *Sterile products are non-returnable.