



## Two Excellent Pediatric Stirrup Choices— E-Z Lift Stirrups & Lithotomy Stirrups



The **Pediatric E-Z Lift Stirrups** are as easy-to-use as its forerunner. Exclusive axial rotation from the handle for superior patient positioning. Boot adjusts with a single clamp which simplifies patient positioning. Boot pads included.

- E-Z Lift feature easily control the patient's leg when moving it into the operative position
- Floating boot automatically adjusts when repositioning the pediatric boot stirrup
- Range: Abduction +25°, Adduction -9°, Lithotomy range of +80° to -30° & externally rotates out 20°

**#800-0342-PUPS** 160 lb (73 kg) patient capacity  
**#800-0342-KIDS** 250 lb (113 kg) patient capacity  
**\$3,644 pr**

**#800-0238** Pediatric E-Z Lift Stirrups Bundle  
 Includes 1 pair of each stirrup **\$6,488 set**

**SpringLoc Clamp:** Required Accessory: Purchase Two  
**#800-0338 US #800-0338-JPN JPN \$199 ea**  
**#800-0338-EU EU #800-0338-SWISS SWISS**  
**#800-0338-DEN DEN**

**Optional Accessories**  
**#800-0074-P** Pedi. E-Z Lift Stirrup Dolly **\$564 ea**  
**#508-1295** PUPS Reusable Boot Pads **\$425 pr**  
**#508-1296** KIDS Reusable Boot Pads **\$475 pr**



The **Pediatric Lithotomy Stirrups** are a comfortable boot stirrup. Specifically designed with the pediatric patient in mind. Pediatric stirrups attach to side rail of any OR table using the Schure Socket XL (sold separately). Boot pads included.

- Pediatric boot stirrup design encapsulates patient's lower leg
- Floating boot automatically adjusts when repositioning the pediatric boot stirrup

**#800-0233** PUPS, 160 lb (73 kg) patient capacity  
**#800-0234** KIDS, 250 lb (113 kg) patient capacity  
**\$1,000 pr**

**#800-0237** Pediatric Lithotomy Stirrups Bundle  
 Includes 1 pair of each stirrup **\$1,600 set**

Required Accessory: Order Two  
**Schure Socket XL**  
**#800-0134 US #800-0134-UK UK**  
**#800-0134-EU EU #800-0134-DEN DEN \$250 ea**

**Optional Accessories**  
**#800-0074-L** Lithotomy Stirrup Dolly **\$564 ea**  
**#508-1295** PUPS Reusable Boot Pads **\$425 pr**  
**#508-1296** KIDS Reusable Boot Pads **\$475 pr**

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351