

## Lithotomy Stirrups

*Effectively positions patients for gynecologic, urologic or laparoscopic procedures*

Patient positioning using the Lithotomy Stirrups features an extended lateral boot design that supports the patient's leg entirely, preventing superficial nerve damage and prevents unwanted leg rotation intra-operatively.

### **Features**

- Use for long or short procedures
- Extended soft boot pads included
- 350 lb. patient weight capacity
- Two required Schure Sockets sold separately
- Sold in pairs

**#800-0049** Lithotomy Stirrups

### **Required Schure Socket: Purchase Two**

- #800-0006** US
- #800-0006-DEN** DEN
- #800-0006-EU** EU
- #800-0006-JPN** JPN
- #800-0006-UK** UK

### **Optional Accessory**

**#800-0074-L** Lithotomy Stirrup Dolly

