



## Extra-Wide Armboard

*Almost twice as wide for a wide range of needs*

This surgical armboard has all the features of our Standard Armboard. SchureMed also offers an adjustable-height, extra-wide armboard that has all the features of the adjustable-height armboard.

### **Features**

- 11" x 26.5"L
- Fits standard U.S. operating room table side rails and is lightweight, extremely durable and radiolucent
- Connects to O.R. table side rails simply and locks securely
- To swing armboard, pull finger-trigger at end of armboard and move to desired position
- Release trigger and it locks tight
- Armboard pads are sold separately

**#800-0018** Extra-Wide Armboard

### **Required Accessory: Choose One**

**#508-0105** 2" (5cm) Deluxe Extra-Wide Pad

**#508-0152** 2" (5cm) Softcare Extra-Wide Pad



**Phone** 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351