



Adj. Height Extra-Wide Armboard

Offers an adjustable height feature

This surgical armboard was created for operating room personnel and accommodates a wide-range of needs. Almost twice as wide, it has all the features of our Standard Armboard. Fits standard U.S. operating room table side rails and is lightweight, extremely durable and radiolucent.



Features

- *Extra-Wide Armboard measures 11" x 26.5"L*
- *Connects to O.R. table side rails simply and locks securely*
- *To swing armboard, pull finger-trigger at end of armboard and move to desired position*
- *Release trigger and it locks tight*

Adjustable Height Extra-Wide Armboard

#800-0017 US

#800-0017-EU EU

#800-0017-UK UK

#800-0017-DEN DEN

Required Accessory: Choose One

#508-0105 2" (5cm) Deluxe Extra-Wide Pad

#508-0152 2" (5cm) Softcare Extra-Wide Pad

Phone 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351