



Schure Slide, Anti-Static, Rigid

Extremely durable patient transfer system

Schure Slide Anti-Static, Rigid is a smooth-surfaced, rigid, and extremely durable patient transfer system. After being placed under a patient, personnel can effortlessly slide the patient from one gurney or table to another without lifting, the prime cause of back injury among hospital staff.

Features

- Reduce back injuries and increase patient comfort
- Non-disposable, no need to waste red bags
- Won't stain – spray and wipe clean
- Speed up turnover time – hang Schure Slide in every O.R.
- Radiolucent and not made with natural rubber latex
- 33"L x 22"W x 3/16" (84cm x 60cm x .5cm)

#800-0068 Schure Slide, Anti-Static, Rigid



Important! Please note that the Schure Slides are made for sliding a patient from one surface to another and are not constructed to support patient weight.