



# SchureMed™ Surgical Stirrup Comparison Chart

	Product	Patient Wt. Limit	Protective Fin Design	Lift Mechanism	Lithotomy Adjustment (From Sterile Field)	Lithotomy Range (From Sterile Field)	Abduction Adjustment (From Sterile Field)	Abduction Range (From Sterile Field)	Lithotomy Indicator
<b>BEST</b>	<b>Great White® HD1000 Stirrups</b> #800-0364	1,000 lbs. (454 kg)	Yes	Yes	Yes All Positions in Range	+75°, -25° Also externally rotates out 20°	Yes All Positions in Range	+25°, -9°	Yes (on rod)
	<b>Great White® Maxima Stirrups</b> #800-0342-M	800 lbs. (363 kg)	Yes	Yes	Yes All Positions in Range	+75°, -30° Also externally rotates out 20°	Yes All Positions in Range	+25°, -9°	Yes (on rod)
	<b>Great White® Platinum Stirrups</b> #800-0342-PL	600 lbs. (272 kg)	Yes	Yes	Yes All Positions in Range	+75°, -30° Also externally rotates out 20°	Yes All Positions in Range	+25°, -9°	Yes (on rod)
	<b>Great White® Robotic Stirrups</b> #800-0342-RF	600 lbs. (272 kg)	Yes	Yes	Yes All Positions in Range	+60°, -55° Also externally rotates out 20°	Yes All Positions in Range	+25°, -9°	Yes (on rod)
<b>BETTER</b>	<b>Great White® Premium Stirrups</b> #800-0342-PR	400 lbs. (181 kg)	No	Yes	Yes All Positions in Range	+75°, -30° Also externally rotates out 20°	Yes All Positions in Range	+25°, -9°	Yes (on rod)
<b>GOOD</b>	<b>Lithotomy Stirrups</b> #800-0049	350 lbs. (159 kg)	Yes	No	No	N/A	N/A	N/A	N/A
	<b>Candy Cane Stirrups</b> #800-0012	350 lbs. (159 kg)	No	No	No	N/A	N/A	N/A	N/A
<b>PEDIATRIC</b>	<b>Pediatric E-Z Lift Stirrups</b> #800-0342-PUPS #800-0342-KIDS	PUPS: 160 lbs. (73 kg) KIDS: 250 lbs. (113 kg)	No	Yes	Yes All Positions in Range	+75°, -30° Also externally rotates out 20°	Yes All Positions in Range	+25°, -9°	Yes (on rod)
	<b>Pediatric Lithotomy Stirrups</b> #800-0233 PUPS #800-0234 KIDS	PUPS: 160 lbs. (73 kg) KIDS: 250 lbs. (113 kg)	No	No	No	N/A	N/A	N/A	N/A

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [schuremed.com](http://schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351