500^{LB}_{PWC}**CamLoc Rail Clamp**

- Fully autoclavable & attaches at any rail location over 1-2 drapes (order Large to use with 3-5 drapes)
- Accommodates both flat and round posts: Flat post dimensions up to 1" x 3/8" (2.5cm x .95cm). Round post dimensions up to 5/8"D (1.6cm).
- Patient weight capacity: 800 lbs (363 kg)

#800-0271

#800-0271-EU


#800-0271-JPN

#800-0271-L

#800-00271-DEN

#800-00271-SWISS

#800-00271-UK

500^{LB}_{PWC}**Deluxe Clamp**

- Accommodates both flat and round posts: Flat post dimensions up to 1" x 3/8" (2.5cm x .95cm). Round post dimensions up to 5/8"D (1.6cm).
- Attaches at rail ends or notches vertically or horizontally
- Patient weight capacity: 720 lbs (327 kg)

#800-0248

#800-0248-EU

#800-0248-SWISS

#800-0248-UK

#800-0248-DEN

#800-0248-JPN

800^{LB}_{PWC}**Schure™ Socket**

- Accommodates round posts up to 5/8" (2cm) in diameter
- Attaches at any rail location
- Patient weight capacity: 800 lbs (363 kg)

#800-0006

#800-0006-EU


#800-0006-JPN

#800-0006-UK

#800-0006-DEN

#800-0006-SWISS

Patented

720^{LB}_{PWC}**Simple Clamp**

- Accommodates flat posts up to 1" x 3/8" (2.5cm x .95cm)
- Attaches at rail ends or notches
- Accessories can be mounted vertically or horizontally
- Patient weight capacity: 720 lbs (327 kg)

#800-0228

#800-0228-EU

#800-0228-SWISS

#800-0228-UK

#800-0228-DEN

#800-0228-JPN

800^{LB}_{PWC}**SpringLoc™ Clamp**

- Accommodates flat posts up to 1" x 1/4" (2.5cm x .6cm)
- Spring-loaded for easy attachment
- Attaches at any rail location
- Patient weight capacity: 800 lbs (363 kg)

#800-0338

#800-0338-EU

#800-0338-JPN

#800-0338-SWISS

350^{LB}_{PWC}**Universal Rail Clamp**

- Accommodates both flat and round posts: Flat post dimensions up to 1" x 3/8" (2.5cm x .95cm). Round post dimensions up to 5/8"D (1.6cm).
- Attaches at any rail location
- Accessories can be mounted vertically or horizontally
- Patient weight capacity: 350 lbs (158 kg)

#800-0085

Patented