



## Extremity Positioner

*Contoured cavity keeps leg or arm in place*

Ideal extremity positioning device for prepping a patient's arm or leg surgical site. Perfect for regular use.

### **Features**

- Great for positioning the lower extremity for foot or ankle procedures
- Sealed seams for easy clean-up
- Dimension: 16"L x 12"W x 7"H

**#508-0353** Extremity Positioner



Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351