



Saphenous Vein Harvesting Pad

Perfect pad to expedite procedure

The Saphenous Vein Harvest Pad holds patient's legs and supports the upper thigh at the optimal position for assessment and harvest of the saphenous vein during coronary bypass surgery.

Features

- *Holds legs externally-rotated and slightly-flexed for minimally-invasive, bridging, or open-harvesting procedures*
- *Pad measures 33"L x 21"W x 7"H*
- *The cover is multi-directional stretch cover*
- *Insides consists of a top layer of conforming foam with a firm, supportive base*
- *Water-resistant, not made with natural rubber latex cover*
- *Sealed, watertight seams for easy clean-up*



#508-0385 Saphenous Vein Harvesting Pad

Phone 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351