



Powered Beach Chair

The only powered surgical beach chair on the market with a self-contained motor

Our double-ball joint head positioning system safely secures the patient's head; even your kyphosis patients and allows patient's head to tilt 30° in every direction. Motion-control technology eliminates manual patient positioning with a simple push-button hand control.

Features

- Keeps the patient and O.R. staff safer and more secure
- Attaches to ANY operating room table with our integrated, adjustable mounting clamps
- 5-Piece Deluxe Pad Set and Lateral Brace included
- 1,350 lbs. of lift securely supports arthroscopy patients up to 500 lbs.
- Shoulder pull-away panels provide excellent surgical exposure
- Articulation range 0-90° in 8 seconds in a slow, controlled movement
- 28"W x 30"D x 47"H, 32 lbs

#800-0004 Powered Beach Chair

Required Accessories: Select One Model

#800-0164 Disp. Face Masks, 12/cs*

#800-0369 Disp. Head & Chin Straps, 12/cs*

**Products are non-returnable*



Optional Accessories

#800-0074 Shoulder Chair Dolly

#800-0091 Lateral Brace

#800-0050 Multi-Axis Arm Positioner

#800-0134 Schure Socket XL



Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351