



Standard Armboard

Offers a full range of motion, swings 180°

Simply pull finger trigger at end of armboard and move to required position. Release trigger and it locks tight. To remove, just pull finger trigger and lift off side rail.

Features

- Simply place armboard on OR table side rail and it snaps locked tight
- Pad sold separately
- 26.5"L x 6"W

Standard Armboard

#800-0008

#800-0008-DEN

#800-0008-EU

#800-0008-JPN

#800-0008-UK

Required Accessory: Choose One

#508-0071 Deluxe Armboard Pad

#508-0151 Softcare Armboard Pad



Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351