



## Adjustable Height Armboard

*Eliminates need for multiple size pads*

Accommodates a wide-range of patients and their procedural needs, as well as adjusting vertically to allow the use of x-ray tops without sacrificing a level surface on the O.R. table.

### Features

- To adjust armboard vertically, simply place armboard on side rail, loosen handle and raise to correct table height and tighten handle
- Simply pull finger-trigger at end of armboard and quickly move to required position. Release trigger to lock tight.
- Pad sold separately
- 6"W x 26.5"L



### Adjustable Height Armboard

#800-0009                      #800-0009-DEN  
#800-0009-EU                #800-0009-JPN  
#800-0009-UK                #800-0009-SWISS

### Required Accessory: Choose One

#508-0071 2" Armboard Pad, Deluxe  
#508-0151 2" Armboard Pad, Softcare

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351