

Schure Foot Lexan

Surgical positioning device for total knee replacement procedures

Easily mounts to operating room side rails using the Schure Socket XL, letting the surgeon replace the entire knee while the patient's foot is held firmly in place. Now available in Radiolucent, clear Lexan®

Features

- Disposable conductive, non-latex foot pad included
- Patient's well leg lies flat
- 12"L x 2"D, 2 lbs
- Schure Socket XL sold separately

#800-0011 Schure Foot, Lexan

Schure Socket XL

#800-0134

#800-0134-EU

#800-0134-JPN

#800-0134-UK

#800-0134-DEN

#800-0134-SWISS

Optional Accessory

#508-0111 Disposable Foot Pads, 12/cs*

*Product is non-returnable.



Phone 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351