



## Leg Prepper

*Stainless steel mounting post adjusts up to 30"H*

Prepare the patient's leg for surgical procedures in a secure and comfortable position with the Leg Prepper from SchureMed.

### **Features**

- *Padded horseshoe cradle holds the patient's leg comfortably during the prepping procedure*
- *Stainless steel 30"H x 5/8"D mounting post rotates, raises, and lowers easily to various prepping positions*
- *Schure Socket XL sold separately*

**#800-0013** Leg Prepper

### **Required Accessory**

**Schure Socket XL**

**#800-0134**

**#800-0134-DEN**

**#800-0134-EU**

**#800-0134-JPN**

**#800-0134-JPN**

**#800-0134-SWISS**

### **Replacement Pad**

**#790-0074** Leg Prepper Pad & Cover



**Phone** 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351