




## End Rest Procedure Table w/Leg

*Adjusts up or down to O.R table level quickly & securely*

SchureMed's End Rest Procedure Table w/Leg, for arm and hand surgery, clamps and locks to standard U.S. operating room tables in seconds.

### Features

- Choose *Lightweight & Radiolucent Carbon Fiber or Phenolic*
- Table measures 16"W (41 cm) at narrow end, 26"L x 34" (66 cm x 86 cm) at wider end
- Radiolucent for all C-arm x-ray images
- Height-adjustable leg w/twist handle locking system
- Deluxe foam pad included
- Two required clamps sold separately
- Replacement Deluxe Foam Pad can be purchased separately



*(Carbon Fiber model shown in photo)*

**#800-0026-ER-P** End Rest Procedure Table w/Leg, Phenolic

**#800-0026-ER-CF** End Rest Procedure Table w/Leg, Carbon Fiber

**Two Clamps Required:** Choose *Universal Clamp* for US or *Deluxe Rail Clamp* for Outside US  
**US - Universal Rail Clamps #800-0085**

**Outside US - Deluxe Rail Clamps**

**#800-0248-EU**

**#800-0248-UK**

**#800-0248-SWISS**

**#800-0248-JPN**

**#800-0248-DEN**

**Optional Accessory**

**#508-0101** End Rest Procedure Table Deluxe Pad

**Phone** 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351