



Hourglass Procedure Table w/Leg

Especially good for mid-arm surgical procedures

Adjusts up or down to be level with the O.R. table quickly and securely. Clamps and locks to standard U.S. operating room tables in seconds.

Features

- Choose *Lightweight & Radiolucent Carbon Fiber or Phenolic*
- *End rail accommodates traction attachment and locks to side rail in seconds*
- *Table measures 15"W (38 cm) at middle x 23" x 32" (58 cm x 81 cm) at table ends*
- *Height-adjustable leg w/twist handle locking system*
- *Deluxe foam pad included*
- *Two required clamps sold separately*
- *Replacement Deluxe Foam Pad can be purchased separately*



(Carbon Fiber model shown in photo)

#800-0026-HG-P Hourglass Procedure Table w/Leg, Phenolic

#800-0026-HG-CF Hourglass Procedure Table w/Leg, Carbon Fiber

Two Clamps Required: Choose *Universal Clamp for US* or *Deluxe Rail Clamp for Outside US*
US - Universal Rail Clamps #800-0085

Outside US - Deluxe Rail Clamps

#800-0248-EU

#800-0248-UK

#800-0248-SWISS

#800-0248-JPN

#800-0248-DEN

Replacement Pad

#508-0106 Hourglass Procedure Table Deluxe Pad

Phone 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351