




Minor Procedure Table

Perfect for all arm and hand procedures

SchureMed's Minor Procedure Table, for arm and hand surgery, clamps and locks to standard U.S. operating room tables in seconds.

Features

- Choose *Lightweight & Radiolucent Carbon Fiber or Phenolic*
- *Adjusts up or down to be level with O.R. table quickly and securely, with no leg attachment*
- *Easy-to-use durable surgical table*
- *Radiolucent for all C-arm x-ray images*
- *Table measures 15"W x 30"L (38 cm x 76 cm)*
- *Two required clamps sold separately*
- *No leg attachment needed*
- *Deluxe foam pad included*



(Carbon Fiber model shown in photo)

#800-0026-MI-P Minor Procedure Table, Phenolic

#800-0026-MI-CF Minor Procedure Table, Carbon Fiber

Two Clamps Required: Choose *Universal Clamp for US* or *Deluxe Rail Clamp for Outside US*
US - Universal Rail Clamps #800-0085

Outside US - Deluxe Rail Clamps

#800-0248-EU

#800-0248-UK

#800-0248-SWISS

#800-0248-JPN

#800-0248-DEN

Replacement Pad

#508-1432 Minor Procedure Table Deluxe Pad

Phone 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351