




Major Procedure Table w/Leg

Extremely durable & lightweight – easy-to-use & radiolucent

SchureMed's Major Procedure Table w/Leg, for arm and hand surgery, clamps and locks to standard U.S. operating room tables in seconds. Adjusts up or down to be level with the O.R. table quickly and securely.

Features

- Choose *Lightweight & Radiolucent Carbon Fiber or Phenolic*
- *Radiolucent for all C-arm x-ray images*
- *Table measures 15"W x 34"L (38 cm x 86 cm)*
- *Height-adjustable leg w/twist handle locking system*
- *1" Deluxe foam pad & Leg included*
- *Two required clamps sold separately*
- *Replacement 1" Deluxe Foam Pad can be purchased separately*



(Carbon Fiber model shown in photo)

#800-0026-MJ-P Major Procedure Table w/Leg, Phenolic

#800-0026-MJ-CF Major Procedure Table w/Leg, Carbon Fiber

Two Clamps Required: Choose *Universal Clamp for US* or *Deluxe Rail Clamp for Outside US*
US - Universal Rail Clamps #800-0085

Outside US - Deluxe Rail Clamps

#800-0248-EU

#800-0248-UK

#800-0248-SWISS

#800-0248-JPN

#800-0248-DEN

Replacement Pad

#508-0073 Major Procedure Table Deluxe Pad

Phone 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351