



Knee Crutches

Extremely soft and comfortable crutch liners prevent superficial nerve damage

Knee Crutches rotate 360° & adjust vertically from 12"-19" (31cm-48cm) to provide precise positioning with an easy-to-use ergonomic locking handle

Features

- Simple to set up and position the patient
- Constructed of the highest quality stainless steel
- Mount crutches simply with Schure Socket XL (sold separately)
- Sold in pairs



#800-0041 Knee Crutches

Schure Socket XL

#800-0134

#800-0134-EU

#800-0134-JPN

#800-0134-UK

#800-0134-DEN

#800-0134-SWISS

Optional Accessory

#508-0154 Knee Crutch Repl. Pads

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351