



Lithotomy Stirrups

Effectively positions patients for gynecologic, urologic or laparoscopic procedures

Patient positioning using the Lithotomy Stirrups features an extended lateral boot design that supports the patient's leg entirely, preventing superficial nerve damage and prevents unwanted leg rotation intra-operatively.

Features

- Use for long or short procedures
- Extended soft boot pads included
- 350 lb. patient weight capacity
- Two required Schure Sockets sold separately
- Sold in pairs

#800-0049 Lithotomy Stirrups

Required Schure Socket: Purchase Two

#800-0006	#800-0006-JPN
#800-0006-DEN	#800-0006-UK
#800-0006-EU	#800-0006-SWISS

Optional Accessory

#800-0074-L Lithotomy Stirrup Dolly

Replacement Pads

#508-1415 Platinum Stirrup Boot Pads, Set



Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351