



Narrow Headrest

Table extension mounts in place of standard headrest

The unique shape improves the surgeon's access at the head of the surgical table during ENT, plastic surgery, neurosurgery or ophthalmic surgical procedures.

Features

- Fits all O.R. tables
- 20"W reduces to 8"W at the head
- Radiolucent headrest includes 2" premium pressure management pad
- Wrist Rest adjusts 14" vertically and inside diameter width adjusts 5" to 11" (not included)
- Two required Simple Clamps sold separately
- Deluxe pad included

#800-0052 Narrow Headrest

#800-0052-EU Narrow Headrest (EU)

Required Simple Clamp: Order Two

#800-0228

#800-0228-EU

#800-0228-UK

#800-0228-SWISS


#800-0228-EU

#800-0228-JPN

Optional Accessories

#800-0086 Wrist Rest


#508-0109 Narrow Headrest Pad



Narrow Headrest



Wrist Rest



Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351