




## Nissen Straps

*Keep your patients comfortable during procedures*

Each leg is wrapped with a 1" thick pad to secure the patient regardless of the angle of the patient's position on the table.

### Features

- *Nissen strap offers the support needed during steep reverse Trendelenburg positioning*
- *20"L x 5"W x 1" thick pad offers complete support*
- *Bariatric model also available at 72"L with 32"L x 7"W x 1" thick pad*
- *Attaches to side rail quickly and securely*
- *Two required Simple Clamps sold separately*
- *Sold in pairs*



**#800-0055** Nissen Straps, 48"L x 2"W, 600 lb. weight capacity

**#800-0118** Bariatric Nissen Straps, 72"L x 2"W, 800 lb. weight capacity

**Required Accessory: Order Two Simple Clamp**

**#800-0228**

**#800-0228-EU**

**#800-0228-JPN**

**#800-0228-UK**

**#800-0228-DEN**

**#800-0228-SWISS**

**Phone** 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351