



SchureFoot, Stainless Steel

Surgical positioning device for total knee replacement procedures

Easily mounts to operating room side rails using the Schure Socket XL, letting the surgeon replace the entire knee while the patient's foot is held firmly in place.

Features

- Disposable conductive, non-latex foot pad included
- Patient's well leg lies flat
- 12"L x 18"H x 2"D (31 cm x 46 cm x 5 cm)
- Schure Socket XL sold separately

#800-0062 SchureFoot, Stainless Steel

Required Accessory

Schure Socket XL

#800-0134

#800-0134-EU

#800-0134-JPN

#800-0134-UK

#800-0134-DEN

#800-0134-SWISS

Optional Accessory

#508-0111 SchureFoot Disposable Pads, 12/cs*

**Product is non-returnable.*



Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351