



## Schure Grip Leg Holder

*Provides versatile & excellent leg positioning*

This surgical leg holder offers optimal patient positioning with dual independent braces that freely rotate 360° and spread from 3" to 11"W to fit any size leg, at any angle.

### Features

- To use, twist one brace handle down, position the patient's leg on the leg holder, twist the handle back up, slide the brace horizontally for a snug fit, and lock in place with convenient twist handle
- Use also for ACL procedures
- Braces measure 12"L x 1.5"D x 8"H
- Mounting post measures 16"H x 5/8"D
- Required Schure Socket XL sold separately
- Stainless steel



**#800-0065** Schure Grip Leg Holder

### Schure Socket XL

**#800-0134**

**#800-0134-EU**

**#800-0134-JPN**

**#800-0134-UK**

**#800-0134-DEN**

**#800-0134-SWISS**

### Optional Accessory

**#508-0103** Disposable Premium Pads, 12/cs\*

\*Product is non-returnable.

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351