



Schure Slide, Anti-Static

Slides easier on the bed sheet or stretcher pad

The Schure Slide, Anti-Static in white enables hospital and operating room personnel to easily slide patient from one bed to another. The Anti-Static version slides even easier on the bed sheet or stretcher pad because of its slippery, high-density, anti-static polyethylene material. The flexible slide bridges the gap between uneven surfaces.

Features

- *Anti-Static, white finish*
- *Non-disposable, no need to waste red bags*
- *Radiolucent and not made with natural rubber latex*
- *Reduce back injuries and increase patient comfort*
- *Speed-up turnover time – hang Schure Slide in every O.R.*
- *Won't stain – spray and wipe clean*
- *33"L x 22"W x 1/16" (84cm x 56cm x .2cm)*
- *Every recovery or emergency room should have one*

#800-0067 Schure Slide, Anti-Static



Important! Please note that the Schure Slides are made for sliding a patient from one surface to another and are not constructed to support patient weight.

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351