



## Schure Slide Full Size, Rigid

*Reduce back injuries & increase patient comfort*

The Full-Size Schure Slide, Rigid enables hospital and operating room personnel to easily slide patient from one bed to another. Each shifter slides easily on the bed sheet or stretcher pad because of its slippery, high-density polyethylene material. The rigid slide bridges the gap between uneven surfaces.

### Features

- Non-disposable, no need to waste red bags
- Radiolucent and not made with natural rubber latex
- Speed-up turnover time – hang Schure Slide in every O.R.
- Unique design has ten hand holes for easy insertion and removal of slide
- Won't stain – spray and wipe clean
- 72"L x 22"W x 3/16" (183cm x 56cm x .5cm)

#800-0069 Schure Slide, Full Size, Rigid

**Important!** Please note that the Schure Slides are made for sliding a patient from one surface to another and are not constructed to support patient weight.

