

Table Width Extensions

Unique design allows any table accessory to be mounted to this extension

The operating room bariatric table width extension mounts on the O.R. table side rail offering the surgeon firm support over the surgical table during upper extremity and thoracic surgical procedures.

Features

- Attaches to any O.R. side rail with mounting Simple Clamp (sold separately)
- Choose from three available sizes
- Raise and lower height of extension to be level with O.R. table; either with or without x-ray top
- Deluxe pad included

#800-0081 Table Width Extension, 11 x 8

#800-0079 Table Width Extension, 15 x 8

#800-0077 Table Width Extension, 20 x 4

Required Simple Clamp: Order One

#800-0228 **#800-0228-EU**

#800-0228-UK **#800-0228-DEN**


#800-0228-JPN **#800-0228-SWISS**

Replacement Pads


#508-0132 Table Width Extension Pad, 11 x 8

#508-0114 Table Width Extension Pad, 15 x 8

#508-0075 Table Width Extension Pad, 20 x 4



11 x 8



15 x 8



20 x 4

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351