



Universal Rail Clamp

Holds all accessories in the vertical or horizontal position

Designed for a wide-range of purposes; affixing most operating room table accessories to most OR table designs and brands.

Features

- Spring-loaded and attaches to OR table side rail anywhere along the rail
- Locks securely with a twist of the ergonomic handle
- Accepts both $\frac{5}{8}$ "D (1.6cm) and 1" x $\frac{3}{8}$ " (2.5cm x .95cm) flat bar mounting posts
- 350 lb (158 kg) patient weight capacity



#800-0085 Universal Rail Clamp

Each clamp sold separately.

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351