

Wrist Tower

Innovative Wrist Tower offers a wide variety of positioning capabilities for wrist procedures

Features

- *Simply dial in your tension...eliminates use of sand bags*
- *You can apply traction by simply dialing the knob from 0-25 lbs (0-11 kgs)*
- *Pronation and supination can be easily controlled by the surgeon's preference and locked in place*
- *Easily attaches to side rail of surgical table with CamLoc Clamp (sold separately)*
- *Improves productivity... eliminates the need for assistant to hold extremity*
- *All components are autoclavable*

#800-0135 Wrist Tower


Required Accessory: Order One CamLoc Clamp*

#800-0271	#800-0271-L
#800-0271-UK	#800-0271-EU
#800-0271-DEN	#800-0271-JPN
#800-0271-SWISS	

**Choose US, Standard Clamp to use with 1-2 drapes or US, Large Clamp for 3-5 drapes*

Optional Accessory Finger Traps, Limited Reuse

Please call customer service for product pricing and specifications.



Phone 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351