



TKR Positioner for Total Knee Replacement

Incredibly innovative knee positioner that offers a completely secure, dynamic-range of patient positions for TKR procedures from SchureMed

This is a huge leap forward for the surgeon and patient alike for total knee replacements. The ball-track system at the heel of the leg holder places the operative leg in flexion, extension, tilt and rotation.

Features

- Track Locking Bracket securely locks the Base Plate to the table with two thumb screws
- Locks in any position with a simple turn of the wrist
- Smooth sliding motion for effortless leg positioning
- Well leg lies flat – eliminating the need for well-leg positioners, sandbags and gel pads
- Easy-to-use and completely rigid once locked
- All components are autoclavable

TKR Positioner for Total Knee Replacement

#800-0141 24" L Track System

#800-0171 30" L Track System

CamLoc Clamp (Required Accessory)

#800-0271 (Use with 1-2 drapes)

#800-0271-L (Use with 3-5 drapes)

#800-0271-JPN


#800-0271-EU

#800-00271-DEN

#800-00271-UK

#800-00271-SWISS

New Improved Design!



Required Accessory

#800-0146

TKR Disp. Sterile TKR Boot Liners, 12/cs*

*Sterile products are non-returnable.

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351