

## Two Excellent Pediatric Stirrup Choices— E-Z Lift Stirrups & Lithotomy Stirrups



The **Pediatric E-Z Lift Stirrups** are as easy-to-use as its forerunner. Exclusive axial rotation from the handle for superior patient positioning. Boot adjusts with a single clamp which simplifies patient positioning. Boot pads included.

- *E-Z Lift feature easily control the patient's leg when moving it into the operative position*
- *Floating boot automatically adjusts when repositioning the pediatric boot stirrup*
- *Range: Abduction +25°, Adduction -9°, Lithotomy range of +75° to -30° & externally rotates out 20°*

**#800-0342-PUPS** 160 lb (73 kg) patient capacity  
**#800-0342-KIDS** 250 lb (113 kg) patient capacity  
**#800-0238** Pediatric E-Z Lift Stirrups Bundle  
*Includes 1 pair of each stirrup*

**SpringLoc Clamp:** *Required Accessory: Purchase Two*  
**#800-0338**      **#800-0338-JPN**  
**#800-0338-EU**   **#800-0338-SWISS**

**Optional Accessories**  
**#800-0074-P** Pediatric Stirrup Dolly  
**#508-1295** PUPS Boot Pads  
**#508-1296** KIDS Boot Pads



The **Pediatric Lithotomy Stirrups** are a comfortable boot stirrup. Specifically designed with the pediatric patient in mind. Pediatric stirrups attach to side rail of any OR table using the Schure Socket XL (sold separately). Boot pads included.

- *Pediatric boot stirrup design encapsulates patient's lower leg*
- *Floating boot automatically adjusts when repositioning the pediatric boot stirrup*

**#800-0233** PUPS, 160 lb (73 kg) patient capacity  
**#800-0234** KIDS, 250 lb (113 kg) patient capacity  
**#800-0237** Pediatric Lithotomy Stirrups Bundle  
*Includes 1 pair of each stirrups*

**Required Accessory: Schure Socket - Order Two**  
**#800-0006**      **#800-0006-SWISS**  
**#800-0006-UK**   **#800-0006-JPN**  
**#800-0006-EU**   **#800-0006-DEN**

**Optional Accessories**  
**#800-0074-L** Lithotomy Stirrup Dolly  
**#508-1295** PUPS Boot Pads  
**#508-1296** KIDS Boot Pads