



Bi-Lateral Leg Prepper

Safely & effectively holds both legs for bi-lateral preps

The unique rotating/pivoting design of the holders allows for unlimited configurations of the lower extremity.

Features

- 500 lb (227 kg) patient weight capacity
- Pads remove for easy cleaning
- Required Schure Socket XL sold separately
- Includes 2" (5cm) Deluxe Pads

#800-0236 Bi-Lateral Leg Prepper

Required Accessory

Schure Socket XL

#800-0134

#800-0134-DEN

#800-0134-EU

#800-0134-JPN

#800-0134-JPN

#800-0134-SWISS

Replacement Pad

#508-0159 Universal Prepper Pad



Phone 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351