



CamLoc Clamp

Attaches to the rail at any location over 1-2 drapes

The CamLoc Clamp is designed for use in the sterile field and fully autoclavable.

Features

- Accepts accessories with round posts up to $\frac{5}{8}$ "D (1.6cm) and flat posts 1" x $\frac{3}{8}$ " (2.5cm x .95cm)
- Attaches along the side rail, without the need of notches
- 500 lb (227 kg) patient weight capacity
- Choose Standard Clamp to use with 1-2 drapes or Large Clamp for 3-5 drapes



CamLoc Clamp

#800-0271

#800-0271-L

#800-0271-EU

#800-0271-JPN

#800-0271-UK

#800-0271-DEN

#800-0271-SWISS

Each clamp sold separately.

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351