

Schure Loc[™] XPS

Your New Surgical Assistant

Schure Loc[™] XPS provides rigid fixation as well as unlimited intraoperative range of motion for shoulder arthroscopy, rotator cuff and total shoulder replacement surgical procedures

Features

- With no foot pedals, hoses, hydraulics or power cords around the operating room table, Schure Loc[™] XPS combines strength and flexibility when intra-operatively adjusting the operative limb
- Single lever releases all six joints for unlimited varied positioning and relocks instantly when released
- Smooth, effortless positioning enhances operative site exposure
- Schure Loc[™] XPS attaches easily to any side rails that are within the following dimensions:
Rail Width Size: .23" - .39" (6 mm - 10 mm)
Rail Height Size: .98" - 1.49" (25 mm - 38 mm)
- Convenient Schure Loc[™] Shoulder Kit includes a sterile drape, foam hand/forearm holder and self-adherent wrap
- 500 lb patient weight capacity

#800-0280 Schure Loc[™] XPS



Optional Accessories

#800-0281 Schure Loc[™] XPS Sterile Kit, 10/cs*
Includes Sterile Forearm Wrap, Sterile Self-Adherent Coban[™] and Sterile Drape (fits whole arm).

*Sterile products are non-returnable.

#790-0039 Schure Loc[™] XPS Adapter

#790-0040 Schure Loc[™] XPS Forearm Attachment

#790-0056 Schure Loc[™] XPS Armboard

#790-0057 Schure Loc[™] XPS Counter Traction Post

#790-0077 Schure Loc[™] XPS TKR Positioner

#790-0075 Schure Loc[™] XPS Elbow/Ankle Positioner

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351