



Adjustable Footboard/Table Extension

New adjustable version allows for a variety of positions

The Adjustable Footboard/Table Extension mounts horizontally as a table extension and vertically as a foot support when using reverse Trendelenburg.

Features

- Includes a 2" Pressure Management Pad
- Simply pull on both locking knobs to lower table extension pad to level of the O.R table pad or pull knobs to raise to a 90° position
- Required SpringLoc Clamps sold separately
- 500 lb. Weight Capacity

#800-0285 Adjustable Footboard/
Table Extension

Required SpringLoc Clamp: Order Two

#800-0338

#800-0338-EU

#800-0338-JPN

#800-0338-SWISS

Optional Accessory

#508-1347 Adjustable Footboard/Table Extension Pad



Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351