



SpringLoc Clamp

Instantly attaches to side rail quickly & smoothly

The newly designed SpringLoc Clamp locks patient positioning devices and surgical equipment with a flat mounting post anywhere along an operating table side rail.

Features

- Securely locks accessories with a flat 1" x ¼" (2.5cm x .6 cm) mounting post
- Attaches easily to the operating table side rail without the need of table notches
- SpringLoc Clamp is easily removed without the need for a release lever
- 1,000 lb. (454 kg) patient weight capacity

SpringLoc Clamp

#800-0338

#800-0338-EU

#800-0338-SWISS

#800-0338-JPN

Each clamp sold separately.



Phone 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351