

Great White[®] Platinum Stirrups

Gynecology, laparoscopy & urology procedures are safer than ever

The Great White[®] Platinum Stirrups provide safe & quick positioning assistance for improved surgical site access while offering exceptionally precise placement.

Features

- Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -30° low to 75° high; Adduction -9° and abduction +25°
- The raised lateral fin boots effectively minimize the possibility of nerve damage
- Included Platinum Boot Pads wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning
- Lithotomy indicators help to accurately reposition legs

#800-0342-PL Great White[®] Platinum Stirrups

Required Accessory: Purchase Two SpringLoc Clamp

#800-0338

#800-0338-EU

#800-0338-JPN

#800-0338-SWISS

Optional Accessories

#800-0074-S Stirrup Dolly

#800-0284 Stirrup Wall Rack

#800-0160 Clear View Leg Drapes, Sterile, 10/cs*

*Sterile products are non-returnable.

Replacement Pads

#508-1415 Platinum Stirrup Boot Pads, Set



600 lb. (272 kg)
Patient Weight
Capacity



Patent
Pending

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351